

# Improvement of Pipetting for a Pooling System:

Authors: H. Rothe<sup>1</sup>, J. Nunes<sup>2</sup>, M. Takeuchi<sup>2</sup>

<sup>1</sup> Apartis Information Management, Bingen, GERMANY  
<sup>2</sup> Roche Molecular Systems, Inc., Pleasanton, CA, USA

*Pooling Wizard* software<sup>3</sup> operating the Microlab® STAR IVD<sup>4</sup> and Microlab® STARlet IVD<sup>4</sup> pipettor as pooling instruments

<sup>3</sup> *Pooling Wizard* is a trademark of Apartis Information Management

<sup>4</sup> Microlab® STAR IVD and/or Microlab® STARlet IVD are registered trademarks of Hamilton Corporation

**apartis**  
Information Management

Head office  
apartis  
Information Management AG  
Gewerbepark Dietersheim  
Im Tiergarten 20  
55411 Bingen  
Germany

US Office  
apartis  
Information Management  
530 Means Street  
Suite 120  
Atlanta, Georgia 30318  
USA


QUALITY WITHOUT COMPROMISE

## Background

Project was to analyze whether pipetting quality can be improved on a pipetting system for pooling, which has pressure measurement implemented (i.e. TADM® - Total Aspiration Dispense Monitoring for the Hamilton Microlab® STAR IVD and Microlab® STARlet IVD).

Used pipetting system already contains pressure control range definitions (Figure 1) in order to detect problematic pipetting events (Figure 2) and mark those pipettings as invalid.

Pressure control requires defining "guard bands", which set the acceptable range within a lower limit curve and a higher limit curve.




## Design

The quality of pressure control depends on the width of the guard bands. Too broad of guard bands may not identify problematic pipettings (Figure 3).

Theoretically a narrow guard band provides higher quality. However, in practice a guard band that is too narrow may invalidate acceptable pipetting events (Figure 4).

The goal has been to identify the optimal settings for each pair of guard bands. The complexity has been given by the number of guard bands needed, i.e. one setting for each different volume, container, aspirate/dispense and tip status (new vs. wet).



The *Pooling Wizard* currently supports 5 different primary pooling methods and 13 different secondary pooling methods for resolution testing (Table 1).


Secondary pool volume is 1ml. Overall this required optimizing 47 different guard bands (Table 2) including the different pressure values for new tips versus wet tips (subsequent pipetting for one primary tube).

## Results

The scope of this project was to evaluate whether pipetting guard bands can be optimized using statistical methods. Therefore multiple pipetting test runs for different scenarios were performed. One possible option was to create a test method for each scenario. The second option was to use the *Pooling Wizard*, which creates all pipetting methods including all single steps on demand. The second option has been chosen by the development team, since this option guarantees 100% similarity of test methods and final methods and so provides higher quality to the end product.

Data was collected in single data bases, one for each pooling run. Each data base contains all curves for a pooling run, i.e. one curve for each step. Data collection happened on different instruments in parallel and at different sites.

One part of the project was to setup a master database, where all single data bases could be uploaded. Statistical evaluation happened afterwards. Apartis developed a *TADM data analyzer* tool, which consists of a central Oracle™ database, a decentral upload tool and a statistical package (Figure 5).



The statistical package has been optimized for curve evaluation and contains various statistical functions for curves: Average curves, minimum curves, maximum curves, moving average curves, standard deviation curves (1x, 2x, 3x) and derivation curves (1st, 2nd).

The master database has been built as a data warehouse with dimensions (e.g. site, instrument, channel, date/time, pooling run) and measures (e.g. pressure curves). The statistical package supports any combination of query parameters. This has been proved as valuable to compare potential differences between sites, instruments, channels, pooling runs and even time impact.


Pool size	PRIMARY POOLING				
	96	48	24	6	1
	with or without Deep Well Plate				
SECONDARY POOLING					
Repeat	96	48	24	6	1
	with or without Deep Well Plate				
2D Pooling	96	48	24	6	1
	with or without Deep Well Plate				
Confirmation	96	48	24	6	1
	with or without Deep Well Plate				
Resolution	96	48	24	6	1
	with or without Deep Well Plate				
Repeat (+) <sup>1)</sup>	96	48	24	6	1
	with or without Deep Well Plate				

Repeat (+) = Repeat positive (checked) test in duplicate

Container	PBT <sup>1)</sup>		SPT <sup>2)</sup>		DEEP WELL PLATE		POOLING PLATE	
	Aspirate	Dispense	Aspirate	Dispense	Aspirate	Dispense	Aspirate	Dispense
Tip	New	Wet	New	Wet	New	Wet	New	Wet
Volume	1000 µl	X	X	X	X	X	X	X
700 µl	-	X	-	-	-	X	X	-
500 µl	X	-	X	X	X	-	-	X
334 µl	X	-	-	-	-	-	-	-
250 µl	X	-	X	X	X	-	-	X
167 µl	X	-	X	X	X	-	-	-
135 µl	X	-	-	-	X	-	-	X
125 µl	X	-	X	X	X	-	-	X
92 µl	X	-	X	X	X	-	-	-

PBT = Primary Blood Tube, SPT = Secondary Pool Tube

## Summary



One sample guard band is used to explain the findings in this publication.

The first iteration uses all measured curves. A reduction of 3% data is performed for each measure point in order to remove spikes. Afterwards Min/Max curves and guard bands with a tolerance of 100 Pa [Pascal] are calculated (Figure 6).

Second iteration calculates the average curve for measure points and uses a tolerance band of 100 Pa again (Figure 7).

Since both iterations show too narrow guard bands during start and end of pipettings, standard deviation (here: 3x STDV) is introduced for third iteration (Figure 8).

Fourth iteration then combines third and second iteration - 3x STDV plus 100 Pa tolerance band (Figure 9). The remaining problem stays with the calculated noise on those guard bands.

Therefore fifth iteration adds a smooth algorithm (Moving average) to calculate the optimized guard bands (Figure 10) for a statistical evaluation approach.

Using guard bands to monitor pipetting pressure is a reliable and well established method to improve pipetting quality. A central database that combines curves from multiple sites, instruments and runs is of essential help. The described statistical method allows for increasingly optimized guard band settings as more data is collected and evaluated.

**Pipetting quality can be improved as a result.**